

**Recombinant Human Annexin A10/ANXA10 Protein**

**Cat.NO.: TP05950**

3th Edition

**Synonyms:**Annexin A10; Annexin-10; Annexin-14; ANXA10; ANX14

**Description:**Annexin A10 (ANXA10) contains four Annexin repeats and is a member of the Annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. It is reported that ANXA10 may be a clinical relevant marker for predicting outcome in both early and advanced stages of bladder cancer.

**Form:**PBS

**Molecular Weight:**37.5 kDa

**Sequences:**Met 1-Tyr324

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.