

Recombinant Human Argininosuccinate Synthase/ASS1 Protein(N-6His)

Cat.NO.: TP05999

3th Edition

Synonyms:Argininosuccinate Synthase; Citrulline--Aspartate Ligase; ASS1; ASS

Description:Argininosuccinate Synthase (ASS1) is an urea cycle enzyme with a tetrameric structure composed of identical subunits. ASS1 is involved in the synthesis of arginine and catalyzes that condensation of citrulline and aspartate to argininosuccinate using ATP. ASS1 is important to the urea cycle as it catalyzes the important second last step in the arginine biosynthetic pathway. ASS1 mainly expressed in periportal hepatocytes, but also in most other body tissues. A deficiency of ASS1 causes citrullinemia (CTLN1), an autosomal recessive disease which is characterized by severe vomiting spells and mental retardation.

Form:PBS

Molecular Weight:42.8 kDa

Sequences:Met 1-Lys412

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.