

Recombinant Human Bone Marrow Proteoglycan/BMPG Protein(C-6His)

Cat.NO.: TP06060

3th Edition

Synonyms:Bone Marrow Proteoglycan; BMPG; Proteoglycan 2; Eosinophil Granule Major Basic Protein; EMBP; MBP; Pregnancy-Associated Major Basic Protein; PRG2; MBP

Description:Bone Marrow Proteoglycan (BMPG) is a secreted protein that contains one C-type lectin domain. BMPG is the predominant constituent of the crystalline core of the eosinophil granule. BMPG is highly expressed in placenta and pregnancy serum. BMPG may be involved in antiparasitic defense mechanisms as a cytotoxin and helminthotoxin. BMPG induces non-cytolytic histamine release from human basophils. In addition, BMPG also participated in antiparasitic defense mechanisms and immune hypersensitivity reactions.

Form:PBS

Molecular Weight:24.6 kDa

Sequences:Leu17-Tyr222

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.