

Recombinant Human Sclerostin/SOST Protein(C-6His)

Cat.NO.: TP06102

3th Edition

Synonyms:Sclerostin;SOST; UNQ2976; PRO7455; PRO7476

Description:Sclerostin, also known as SOST, is a member of the Cerberus/DAN family of BMP antagonists. SOST is a secreted glycoprotein with a C-terminal cysteine knot-like (CTCK) domain. It shows sequence similarity to the DAN (differential screening-selected gene aberrative in neuroblastoma) family of bone morphogenetic protein (BMP) antagonists. SOST is produced by the osteocyte and has anti-anabolic effects on bone formation. It is a negative regulator of bone growth that acts through inhibition of Wnt signaling and bone formation.

Form:PBS

Molecular Weight:22.3 kDa

Sequences:Gln24-Tyr213

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.