

Recombinant Human Carbonic Anhydrase 4/CA4 Protein(C-6His)

Cat.NO.: TP06113

3th Edition

Synonyms:Carbonic Anhydrase 4; Carbonate Dehydratase IV; Carbonic Anhydrase IV; CA-IV; CA4

Description:Carbonic Anhydrase 4 (CA4) belongs to the alpha-carbonic anhydrase family. Alpha-carbonic anhydrase is a large family of zinc metalloenzymes that catalyze the reversible hydration of carbon dioxide. Carbonic anhydrase 4 is a glycosylphosphatidyl-inositol-anchored membrane isozyme expressed on the luminal surfaces of pulmonary (and certain other) capillaries and proximal renal tubules. Carbonic anhydrase 4 may stimulate the sodium/bicarbonate transporter activity of SLC4A4 that acts in pH homeostasis. It may have a role in inherited renal abnormalities of bicarbonate transport. Furthermore, Carbonic anhydrase 4 is essential for acid overload removal from the retina and retina epithelium and acid release in the choriocapillaris.

Form:PBS

Molecular Weight:31.4 kDa

Sequences:Ala19-Lys283

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.