

**Recombinant Human MANSC Domain-Containing Protein 1/MANSC1 Protein(C-6His)**

**Cat.NO.: TP06119**

3th Edition

**Synonyms:**MANSC Domain-Containing Protein 1; Loss of Heterozygosity 12 Chromosomal Region 3 Protein; MANSC1; LOH12CR3

**Description:**MANSC domain-containing protein 1(MANSC1) is encoded by MANSC1 gene. MANSC1 is a Single-pass type I membrane protein which contains 1 MANSC domain. It is widely expressed in many tissues and mainly located in membrane.

**Form:**PBS

**Molecular Weight:**39.8 kDa

**Sequences:**Gln27-Leu385

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.