



Recombinant Human Serpin B6/Placental Thrombin Inhibitor Protein(N-Trx, 6His)

Cat.NO.: TP06151

3th Edition

Synonyms:Serpin B6; Cytoplasmic Antiproteinase; CAP; Peptidase Inhibitor 6; PI-6; Placental Thrombin Inhibitor; SERPINB6; PI6; PTI

Description:Serpin B6 belongs to the serpin family. Serpin B6 localizes to the cytoplasm. Serpin B6 is expressed in many tissues, abundantly by mast cells in different tissues and mastocytoma lesions. Serpin B6 may be involved in the regulation of serine proteinases present in the brain or extravasated from the blood. In addition, Serpin B6 may play an important role in the inner ear in the protection against leakage of lysosomal content during stress and loss of this protection results in cell death and sensorineural hearing loss.

Form:PBS

Molecular Weight:60.7 kDa

Sequences:Met 1-Pro376

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.