

Recombinant Human Myozenin-2/MYOZ2 Protein(C-6His)

Cat.NO.: TP06186

3th Edition

Synonyms:Myozenin-2; Calsarcin-1; FATZ-Related Protein 2; MYOZ2; C4orf5

Description:Myozenin 2 (MYOZ2) is a 264 amino acid protein that belongs to the myozenin family. MYOZ2 binds to Calcineurin, a phosphatase that is involved in calcium-dependent signal transduction in diverse cell types. MYOZ2 is one of the sarcomeric proteins and plays an important role in myofibrillogenesis and the modulation of calcineurin signaling. It may serve as intracellular binding proteins involved in linking Z line proteins such as alpha-actinin, gamma-filamin, TCAP/telethonin, LDB3/ZASP and plays an important role in the modulation of calcineurin signaling. Defects in MYOZ2 are the cause of familial hypertrophic cardiomyopathy type 16 (CMH16), a hereditary heart disorder.

Form:PBS

Molecular Weight:30.9 kDa

Sequences:Met 1-Leu264

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.