

Recombinant Mouse Complement factor D Protein

Cat.NO.: TP06198

3th Edition

Synonyms: Complement factor D; 28 kDa adipocyte protein; Adipsin; C3 convertase activator; Properdin factor D; Cfd; Adn; Df

Description: Complement factor D, also known as adipsin, is a member of the chymotrypsin family of serine proteases, which plays an essential role in host defense as the rate-limiting enzyme in the alternative pathway of complement activation. Complement factor D activates a convertase (C3bBb) responsible for cleavage of the complement protein C3, which leads to the activation of terminal complement component C5-9 to form the membrane attack complex on microbial or cellular surfaces. It also functions in the regulation of systemic energy balance and physiologic and pathologic processes, including immunity and inflammation.

Form: PBS

Molecular Weight: 26.5 kDa

Sequences: Ile26-Ser259

Purity: > 95% by HPLC

Concentration:

Endotoxin Level: <1.0 EU per 1 ug of protein (determined by LAL method)

Storage: Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.