

Recombinant Mouse C-C motif Chemokine 8/CCL8/MCP-2 Protein(C-6His)

Cat.NO.: TP06200

3th Edition

Synonyms:C-C motif chemokine 8;Ccl8;Monocyte chemoattractant protein 2;Monocyte chemotactic protein 2;MCP-2;Small-inducible cytokine A8;Mcp2; Scya8

Description:Chemokine ligand 8 (CCL8,MCP-2), is a small secreted cytokine which belongs to the intercrine beta (chemokine CC) family. CCL8 Chemotactic factor attracts monocytes. It can bind heparin.CCL8 functions to activate different immune cells, including mast cells, eosinophils and basophils which are involved in allergic responses, monocytes, and T cells and NK cells which are involved in the inflammatory response. Its ability achieves by binding to different cell surface receptors termed chemokine receptors including CCR1, CCR2B and CCR5. It has been reported that CCL8 is a potent inhibitor of HIV-1 by virtue of its binding to CCR5 which is one of the major co-receptors for HIV-1.

Form:PBS

Molecular Weight:9.8 kDa

Sequences:Glu20-Pro97

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.