

Recombinant Human Resistin Protein (Fc Tag)

Cat.NO.: TP06217

3th Edition

Synonyms:ADSF;FIZZ3;RETN1;RSTN;XCP1

Description: Resistin is an adipocytokine, which has been studied for its role in insulin resistance and recently in inflammation. The RETN and CAP1 polymorphisms and gene expression may be potential biomarkers for breast cancer risk. Resistin (RETN), recently found to be relevant to inflammation and inflammatory disorders.

Form:PBS

Molecular Weight:38.2 kDa

Sequences:Ser17-Pro108

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.