

**Recombinant Human RAB27B Protein (His tag)**

**Cat.NO.: TP06282**

3th Edition

**Synonyms:**C25KG;RAB27B

**Description:**Members of the Rab protein family, including RAB27B, are prenylated, membrane-bound proteins involved in vesicular fusion and trafficking (Chen et al., 1997 [PubMed 9066979]).

**Form:**PBS

**Molecular Weight:**26.6 kDa

**Sequences:**Thr 2-Lys215

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.