

Recombinant Human RAB27B Protein (Fc tag)

Cat.NO.: TP06283

3th Edition

Synonyms:C25KG;RAB27B

Description:Members of the Rab protein family, including RAB27B, are prenylated, membrane-bound proteins involved in vesicular fusion and trafficking (Chen et al., 1997 [PubMed 9066979]).

Form:PBS

Molecular Weight:50.8 kDa

Sequences:Thr 2-Lys215

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.