

Recombinant Human TWF1 / PTK9 / Twinfilin-1 Protein (His & GST tag)

Cat.NO.: TP06438

3th Edition

Synonyms:A6;MGC23788;MGC41876;PTK9

Description:Twinfilin-1, also known as Protein A6, Protein tyrosine kinase 9, TWF1 and PTK9, is a cytoplasm protein which belongs to the actin-binding proteins ADF family and Twinfilin subfamily. Twinfilin-1 (TWF1 / PTK9) is a highly conserved actin monomer-binding protein that regulates cytoskeletal dynamics in organisms from yeast to mammals. In addition to the mammalian twinfilin-1, a second protein with approximately 65% sequence identity to twinfilin-1 exists in mouse and humans. TWF1 / PTK9 is expressed at high levels in the colon, testis, ovary, prostate and lung. It is expressed at lower levels in the brain, bladder and heart. It is not detected in liver. TWF1 / PTK9 is an actin-binding protein involved in motile and morphological processes. It inhibits actin polymerization, likely by sequestering G-actin. By capping the barbed ends of filaments, it also regulates motility. TWF1 / PTK9 seems to play an important role in clathrin-mediated endocytosis and distribution of endocytic organelles.

Form:PBS

Molecular Weight:57 kDa

Sequences:Met 1-Asp 252

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.