

Recombinant Mouse PGRL / IGSF8 Protein (Fc tag)

Cat.NO.: TP06476

3th Edition

Synonyms:AA033172;AU041109;ESTM34;EWI-2;KCT-4;PGRL;R75384

Description:

Form:PBS

Molecular Weight:85.8 kDa

Sequences:Met 1-Thr 577

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.