

**Recombinant Human PTPN12 Protein (aa 1-355, His & GST tag)**

**Cat.NO.: TP06489**

3th Edition

**Synonyms:**PTP-PEST;PTPG1

**Description:**PTPN12 is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. PTPN12 contains a C-terminal PEST motif, which serves as a protein–protein interaction domain, and may be related to protein intracellular half-life. PTPN12 was found to bind and dephosphorylate the product of oncogene c-ABL, thus may play a role in oncogenesis. PTPN12 was shown to interact with, and dephosphorylate, various of cytoskeleton and cell adhesion molecules, such as p130 (Cas), CAKbeta/PTK2B, PSTPIP1, and paxillin, which suggested its regulatory roles in controlling cell shape and mobil.

**Form:**PBS

**Molecular Weight:**69.4 kDa

**Sequences:**Met 1-Gln355

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.