

Recombinant Mouse CHODL / Chondrolectin Protein (His tag)

Cat.NO.: TP06589

3th Edition

Synonyms:3110074E07Rik;MT75;PRED12

Description:

Form:PBS

Molecular Weight:23.4 kDa

Sequences:Met 1-Asn 216

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.