

**Recombinant Mouse CHODL / Chondrolectin Protein (Fc tag)**

**Cat.NO.: TP06590**

3th Edition

**Synonyms:**3110074E07Rik;MT75;PRED12

**Description:**

**Form:**PBS

**Molecular Weight:**48.9 kDa

**Sequences:**Met 1-Asn 216

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.