



Recombinant Mouse SerpinB12 Protein (His tag)

Cat.NO.: TP06617

3th Edition

Synonyms:2300003F07Rik;4833409F13Rik

Description:Serpins are the largest and most diverse family of serine protease inhibitors which are involved in a number of fundamental biological processes such as blood coagulation, complement activation, fibrinolysis, angiogenesis, inflammation and tumor suppression and are expressed in a cell-specific manner. Most serpins are secreted and attain physiologic concentrations in the blood and extracellular fluids. Mouse SerpinB12 is a cytoplasm protein which belongs to the serpin family and Ov-serpin subfamily. It is expressed in many tissues, including brain, bone marrow, lymph node, heart, lung, liver, pancreas, testis, ovary, and intestine. SerpinB12 inhibits trypsin and plasmin, but not thrombin, coagulation factor Xa, or urokinase-type plasminogen activator.

Form:PBS

Molecular Weight:49.3 kDa

Sequences:Met 1-Pro 423

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.