



Recombinant Human ITGA8 & ITGB1 Heterodimer Protein

Cat.NO.: TP06690

3th Edition

Synonyms:ITGA8;RP11-479G22.2;CD29;FNRRB;GPIIA;MDF2;MSK12;VLA-BETA;VLAB

Description:

Form:PBS

Molecular Weight:196 (113 + 83) kDa

Sequences:Met 1-Leu 1012

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.