

Recombinant Human Carbonic Anhydrase IV / CA4 Protein (His tag)

Cat.NO.: TP06692

3th Edition

Synonyms:CAIV;Car4;RP17

Description:The carbonic anhydrases (or carbonate dehydratases) are classified as metalloenzyme for its zinc ion prosthetic group and form a family of enzymes that catalyze the rapid interconversion of carbon dioxide and water to bicarbonate and protons, a reversible reaction that takes part in maintaining acid-base balance in blood and other tissues. The carbonic anhydrase I (CA) family consists of at least 11 enzymatically active members and a few inactive homologous proteins. Carbonic anhydrase IV (CAIV) is a membrane-associated enzyme anchored to plasma membrane surfaces by a phosphatidylinositol glycan linkage. CAIV is a high-activity isozyme in CO₂ hydration comparable to that of CAII. Furthermore, CAIV is more active in HCO₃⁻ dehydration than is CAII. However, the esterase activity of CAIV is decreased 150-fold compared to CAII.

Form:PBS

Molecular Weight:31.7 kDa

Sequences:Met 1-Lys283

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.