



Recombinant Human ITGA6 & ITGB1 Heterodimer Protein

Cat.NO.: TP06694

3th Edition

Synonyms:CD49f;ITGA6B;VLA-6

Description:

Form:PBS

Molecular Weight:198 (115 + 83) kDa

Sequences:Met 1-Ser 1011

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.