

Recombinant Human KIAA1279 Protein (His & GST tag)

Cat.NO.: TP06765

3th Edition

Synonyms:KBP;KIAA1279;TTC20

Description:

Form:PBS

Molecular Weight:99.6 kDa

Sequences:Met 1-Thr621

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.