

Recombinant Mouse COLEC10 Protein (Fc tag)

Cat.NO.: TP06775

3th Edition

Synonyms:CL-L1

Description:COLEC10 is a member of the C-lectin family. Members of this family possess collagen-like sequences and carbohydrate recognition domains. The other members of this family are secreted proteins and bind to carbohydrate antigens on microorganisms facilitating their recognition and removal. This gene product is a cytosolic protein, a characteristic that suggests that it may have different biological functions than other C-lectins.

Form:PBS

Molecular Weight:47 kDa

Sequences:Cys119-Lys277

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.