

Recombinant Human RAB11B Protein (His Tag)

Cat.NO.: TP06807

3th Edition

Synonyms:H-YPT3

Description:The two recurrent dominant mutations in RAB11B leading to a neurodevelopmental syndrome, likely caused by altered GDP/GTP binding that inactivate the protein and induce GEF binding and protein mislocalization.

Form:PBS

Molecular Weight:26.3 kDa

Sequences:Gly2-Cys215

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.