



Recombinant Human RAB31 / Rab22B Protein (His tag)

Cat.NO.: TP06818

3th Edition

Synonyms:Rab22B

Description:Small GTP-binding proteins of the RAB family, such as RAB31, play essential roles in vesicle and granule targeting (Bao et al., 2002 [PubMed 11784320])

Form:PBS

Molecular Weight:23.9 kDa

Sequences:Met 1-Cys195

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.