

Recombinant Human GMPR Protein (His Tag)

Cat.NO.: TP06827

3th Edition

Synonyms:GMPR1

Description:GMPR, also known as GMPR1, belongs to the IMPDH/GMPR family. This family of enzymes includes IMP dehydrogenase and GMP reductase. These enzymes are involved in purine metabolism and adopt a TIM barrel structure. GMPR is an enzyme that catalyzes the irreversible and NADPH-dependent reductive deamination of GMP to IMP. GMPR functions in the conversion of nucleobase, nucleoside and nucleotide derivatives of G to A nucleotides, and in maintaining the intracellular balance of A and G nucleotides.

Form:PBS

Molecular Weight:39.2 kDa

Sequences:Met 1-Ser345

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.