

Recombinant Mouse GLIPR1 / RTVP1 Protein (His tag)

Cat.NO.: TP06860

3th Edition

Synonyms:2410114O14Rik;mRTVP-1;RTVP-1;RTVP1

Description:Glioma pathogenesis-related protein 1, also known as Protein RTVP-1, GLIPR1 and GLIPR, is a single-pass membrane protein which belongs to the CRISP family. GLIPR1 / RTVP-1 was expressed in high levels in glioblastomas, whereas its expression in low-grade astrocytomas and normal brains was very low. Transfection of glioma cells with small interfering RNAs targeting GLIPR1 / RTVP-1 decreased cell proliferation in all the cell lines examined and induced cell apoptosis in some of them. Overexpression of GLIPR1 / RTVP-1 increased astrocyte and glioma cell proliferation and the anchorage-independent growth of the cells. In addition, overexpression of GLIPR1 / RTVP-1 rendered glioma cells more resistant to the apoptotic effect of tumor necrosis factor-related apoptosis-inducing ligand and serum deprivation. GLIPR1 / RTVP-1 regulated the invasion of glioma cells was evident by their enhanced migration through Matrigel and by their increased invasion in a spheroid confrontation assay. The increased invasive potential of the GLIPR1 / RTVP-1 overexpressors was also shown by the increased activity of matrix metalloproteinase 2 in these cells. The expression of GLIPR1 / RTVP-1 is correlated with the degree of malignancy of astrocytic tumors and that GLIPR1 / RTVP-1 is involved in the regulation of the growth, survival, and invasion of glioma cells. GLIPR1 / RTVP-1 is a potential therapeutic target in gliomas.

Form:PBS

Molecular Weight:25.1 kDa

Sequences:Met1-Thr223

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.