

Recombinant Human Troponin C / TNNC1 Protein

Cat.NO.: TP06907

3th Edition

Synonyms:CMD1Z;CMH13;TN-C;TNC;TNNC

Description:Troponin I, also known as TNNC1, is part of the troponin complex. This complex contains 3 subunits: troponin I (TnI), troponin T (TnT) and troponin C (TnC). Troponin I is the inhibitory subunit, blocking actin-myosin interactions and thereby mediating striated muscle relaxation. It binds to actin in thin myofilaments to hold the actin-tropomyosin complex in place. Because of it myosin cannot bind actin in relaxed muscle. When calcium binds to the Troponin C it causes conformational changes which lead to dislocation of troponin I and finally tropomyosin leaves the binding site for myosin on actin leading to contraction of muscle.

Form:PBS

Molecular Weight:18.4 kDa

Sequences:Met 1-Glu161

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.