

**Recombinant Human Troponin C / TNNC1 Protein**

**Cat.NO.: TP06907**

3th Edition

**Synonyms:**CMD1Z;CMH13;TN-C;TNC;TNNC

**Description:**Troponin I, also known as TNNC1, is part of the troponin complex. This complex contains 3 subunits: troponin I (TnI), troponin T (TnT) and troponin C (TnC). Troponin I is the inhibitory subunit, blocking actin-myosin interactions and thereby mediating striated muscle relaxation. It binds to actin in thin myofilaments to hold the actin-tropomyosin complex in place. Because of it myosin cannot bind actin in relaxed muscle. When calcium binds to the Troponin C it causes conformational changes which lead to dislocation of troponin I and finally tropomyosin leaves the binding site for myosin on actin leading to contraction of muscle.

**Form:**PBS

**Molecular Weight:**18.4 kDa

**Sequences:**Met 1-Glu161

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.