

Recombinant Human TRP1 / TYRP1 Protein (His tag)

Cat.NO.: TP06910

3th Edition

Synonyms:b-PROTEIN;CAS2;CATB;GP75;OCA3;TRP;TRP1;TYRP

Description:Tyrosinase-related protein 1, also known as TYRP1 or TRP1, is a melanosomal enzyme that belongs to the tyrosinase family and plays an important role in the melanin biosynthetic pathway. Mutations in this enzyme are the cause of rufous oculocutaneous albinism and oculocutaneous albinism type III. TYRP1 / TRP1 is involved in the oxidation of 5,6-dihydroxyindole-2-carboxylic acid (DHICA) into indole-5,6-quinone-2-carboxylic acid. This enzyme may regulate or influence the type of melanin synthesized. The expression of Tyrosinase-related protein 1 (TYRP1) is regulated by the microphthalmia-associated transcription factor (MITF). There is mounting evidence demonstrating that in addition to its role in eumelanin synthesis, TYRP1 is involved in maintaining stability of tyrosinase proliferation and melanocyte cell death.

Form:PBS

Molecular Weight:52 kDa

Sequences:Met 1- Arg 471

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.