

**Recombinant Mouse NKR-P1A / Klrb1a Protein (His tag)**

**Cat.NO.: TP06914**

3th Edition

**Synonyms:**ly-55A;Ly55a;NKR-P1.7;NKR-P1A;Nkrp1-a;NKRP12;Nkrp1a

**Description:**

**Form:**PBS

**Molecular Weight:**21 kDa

**Sequences:**Gln 67-His 227

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.