

**Recombinant Mouse CD79A / Ig-alpha Protein (ECD, Fc tag)**

**Cat.NO.: TP06977**

3th Edition

**Synonyms:**lg-alpha;lga;lgalpha;Ly-54;Ly54;mb-1

**Description:**

**Form:**PBS

**Molecular Weight:**39 kDa

**Sequences:**Met1-Arg137

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.