

**Recombinant Rat CLEC4A2 / DCIR Protein (His tag)**

**Cat.NO.: TP06985**

3th Edition

**Synonyms:**CLEC4A2;Dcir2

**Description:**

**Form:**PBS

**Molecular Weight:**21.8 kDa

**Sequences:**Gln70-Ser235

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.