

**Recombinant Human LRG1 Protein (Fc tag)**

**Cat.NO.: TP07022**

3th Edition

**Synonyms:**HMFT1766;LRG

**Description:**LRG1 belongs to the leucine-rich repeat (LRR) family. Members of this family are involved in protein-protein interaction, signal transduction, and cell adhesion and development. LRG1 is expressed during granulocyte differentiation. It contains 4 LIM zinc-binding domains and 1 Rho-GAP domain.

**Form:**PBS

**Molecular Weight:**61.3 kDa

**Sequences:**Met 1-Gln347

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.