

**Recombinant Human IL-12 (IL12A & IL12B Heterodimer) Protein (His Tag)**

**Cat.NO.: TP07032**

3th Edition

**Synonyms:**CLMF;CLMF2;IL-12A;IL-12B;IL12;IL12 p70;IMD28;IMD29;NFSK;NKSF;NKSF1;NKSF2;P35

**Description:**

**Form:**PBS

**Molecular Weight:**59.6 kDa

**Sequences:**Met 1-Ser328

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.