

Recombinant Human MKK6 / MEK6 / MAP2K6 Protein (207 Ser/Asp, 211 Thr/Asp, His & GST tag)

Cat.NO.: TP08434

3th Edition

Synonyms:MAPKK6;MEK6;MKK6;PRKMK6;SAPKK-3;SAPKK3

Description:Dual specificity mitogen-activated protein kinase kinase 6, also known as MAP kinase kinase 6, MAPKK 6, MAPK / ERK kinase 6, SAPKK3, MAP2K6 and MKK6, is a protein which belongs to the protein kinase superfamily, STE Ser / Thr protein kinase family and MAP kinase kinase subfamily. MAP2K6 / MKK6 contains one protein kinase domain. Mitogen-activated protein kinases are members of a conserved cascade of kinases involved in many signal transduction pathways. They stimulate phosphorylation of transcription factors in response to extracellular signals such as growth factors, cytokines, ultraviolet light, and stress-inducing agents. MAP2K6 / MKK6 exists in a variety of alternatively spliced isoforms with distinct patterns of tissue expression. Isoform 2 of MAP2K6 / MKK6 is only expressed in skeletal muscle. Isoform 1 of MAP2K6 / MKK6 is expressed in skeletal muscle, heart, and in lesser extent in liver or pancreas.

Form:PBS

Molecular Weight:65.3 kDa

Sequences:Met 1-Asp 334, Ser 207/Asp, Thr 211/Asp

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.